

# Biomedical Applications of Molecular Spectroscopy

Mike Kayat


B&W Tek, Inc  
19 Shea Way  
Newark, DE 19713  
United States of America  
+1 302 368 7824

[mikek@bwtek.com](mailto:mikek@bwtek.com)

## Overview

- Molecular spectroscopy is a large and expanding discipline.
- One of the key areas that is driving this expansion is biomedical applications of molecular spectroscopy (“biospectroscopy”)
- Low cost compact modular spectrometers are enabling more applications
- Several methods are used in the analysis of DNA, proteins, viruses, bacteria, cells, tissues and these mainly use ultra-violet (UV), visible, near infrared (NIR) spectrometers together with a range of sampling techniques
- A review of biomedical applications using these different spectrometers, particularly for UV absorption, NIR transmission and reflectance, together with Raman scattering, will be discussed
- Rapid adoption of spectrometers which are now being deployed for biomedical analysis, for example: antioxidants, arteries particularly atherosclerosis, bone, cancer tumors, diabetes, fertility, infections and other examples

## Biospectroscopy – Large Area of Molecular Spectroscopy




Most biomedical applications based on:

- UV absorption spectroscopy
- NIR spectroscopy
- Raman spectroscopy

Source: SDI 2011


## Biospectroscopy – Large Area of Molecular Spectroscopy




- Biological macromolecules:
  - nucleic acids, proteins, lipids
- Blood disorders:
  - anemias, leukemias, thalassemias
- Cancer diagnosis:
  - brain, breast, cervical, colon and other
- Chemical processes in live blood cells:
  - malaria, drug reactions
- DNA in chromosomes, pigment in granulocytes, RBCs, hepatocytes
- Immuno assays
- Organelles, cells, micro-organisms, bacteria, phytoplankton, neurotoxins, viruses
- Tissue analysis:
  - Alzheimer's, artery, breast, bone, cervix, embryo media, esophagus, gastro-intestinal tract, prostate

# Biospectroscopy – Biophysics

Photo-Molecular Interactions



## Biospectroscopy – Examples of UV Absorption Spectroscopy - $\mu\text{L}$



For nucleic acid quantification, the Beer-Lambert equation is:

$$c = (A * e)/b$$

Where  $c$  is the nucleic acid concentration in  $\text{ng}/\mu\text{L}$ ,  $A$  is the absorbance in AU,  $e$  is the wavelength-dependent extinction coefficient in  $\text{ng}\cdot\text{cm}/\mu\text{L}$  and  $b$  is the path length in cm. The generally accepted extinction coefficients for nucleic acids are:

- Double-stranded DNA: 50
- Single-stranded DNA: 33
- RNA: 40

Source: Nanodrop Technologies

- Acquire UV-visible absorption spectrum between 220 nm to 750 nm
- Measure peak ratios 260 nm/280 nm; for DNA, the peak should be at 260 nm, the 260/280 ratio should be between 1.8 and 2.0
- Concentration measurements of DNA, RNA, dyes, proteins, cell cultures
- Typically 1 pg per  $\mu\text{L}$  sensitivity
- Each sample is measured using two different path lengths (1mm and 0.2 mm), providing a wide dynamic range ( $\sim 2 \text{ ng}/\mu\text{L}$  - 3700  $\text{ng}/\mu\text{L}$  dsDNA)
- The short path length provides 50 times higher in concentration than can be measured on classical 1 cm cuvette-based systems

## Biospectroscopy – Advantages of Vibrational Spectroscopy

- Vibrational spectroscopy including Raman scattering and near infrared(NIR) absorption, transmission and reflectance has large and increasing potential
- Number of degrees of freedom in a complex molecule is  $3N - 6$ , where  $N$  is the number of atoms and for a large molecule like a biopolymer with  $\sim O(100)$  atoms, there are effectively three degrees of freedom for each atom
- Spectra tend to be complex with many lines, but different types of molecules have characteristic behavior
- Different distinct spectra for tissues and pathologies
- Objective is to identify spectral characteristics that are associated with a various states including disease


## Biospectroscopy – Comparison Raman vs IR

	Near-IR	Mid-IR	Raman
Spectral range (cm <sup>-1</sup> )	13,300–3300	4000–400	4000–50
Analysis of:			
Gases	No	Yes	Yes
Liquids	Yes	Yes	Yes
Solids	Yes	Yes	Yes
Aqueous systems	Difficult	Very difficult	Yes
Macroscopic samples	Yes	Yes	Yes
Microscopic samples	No	Yes	Yes
Signal	Strong	Strong	Weak
Sampling	Easy	Difficult	Easy
Through glass windows	Yes	No	Yes
In situ	No	No	Yes
Quantitative	Yes	Difficult	Yes
Noninvasive	Yes	No	Yes
Fiber optic interfacing	Yes	No	Yes
Information content	Low. Limited to O–H, N–H, and C–H vibrations	High	High
Reaction monitoring and modeling	Requires chemometrics	Yes	Yes



# Biospectroscopy – Examples of NIR Spectroscopy – IVF

- In vitro fertilization application
- Embryo viability testing



NIR-based spectrometer with touch screen  
25 cm (H) x 28 cm (W) x 42 cm (D)

Three components:

1. Analysis instrument
2. Sample cells, dark cell & reference cells
3. Temperature stabilizer (23.8°C)




Single use sample cell with RFID  
5.4 cm (H) x 2.2 cm (W) x 0.6 cm (D)

Source: Molecular Biometrics

## Biospectroscopy – Examples of NIR Spectroscopy - IVF

- Metabolomic profile of spent embryo culture (2 -3 days)
- Target specific Oxidation Stress (OS) biomarkers in embryonic culture media and measure NIR absorption
- Identify metabolomic differences in viable compared with non-viable embryos
- Non-invasive spectroscopic method for prediction of implantation potential of embryos in IVF
- Measure changes in the -CH, -NH, -CH, -SH functional groups between 900-1700 nm
- **Viability Score** indicates the SET embryo reproductive potential
- Established correlation between Viability Score and fetal cardiac activity within 12 weeks



Source: Molecular Biometrics


## Biospectroscopy – Advantages of Raman Spectroscopy

- Wavelength selection
- No water interference (physiological state)
- No biopsy required
- Directly measures molecules
- Small concentrations
- Chemical composition
- Morphological analysis
- Quantitative analysis from sharp spectral peaks
- *In vivo* diagnosis
- High spatial resolution
- Raman can distinguish numerous pathologies
- Raman probes can distinguish disease
- Raman can now be used in deeper tissues

## Biospectroscopy – Examples of Raman spectroscopy

- Raman spectroscopy ‘fingerprints’ molecules by characterizing interactions between photons and molecular vibrations (*unique for each biomolecule*)
- Near-infrared excitation is preferred for biomedical applications
- Recent optical fiber probe developments allow accurate real-time analysis *in vivo*
- Alternative to immunofluorescence staining (*clinical diagnostics*)
- New areas of research are promising for widespread clinical applications (*cell biology, imaging, tissue engineering, pathogens, pharmacology*)


# Biospectroscopy – Examples of Raman Spectroscopy




Hanlon *et al.* "Prospects for *in vivo* Raman spectroscopy," *Phys Med Biol* **45**: R1 (2000)

# Biospectroscopy – Examples of Raman Spectroscopy

## Raman Spectral Pathology of Atherosclerosis




## Biospectroscopy – Examples of Raman Spectroscopy



Buschman HPJ, et al. *Cardiovascular Pathology* 10(2), 69-82 (2001)

# Biospectroscopy – Examples of Raman Spectroscopy




<sup>1</sup>Gloucestershire Royal Hospital, Gloucester, UK

100 mW excitation, 785 nm, 1 second collection





# Biospectroscopy – Examples of Raman spectroscopy





# Biospectroscopy – Examples of Raman spectroscopy

Application	Raman spectroscopy	Market potential
Biological macromolecules (nucleic acids, proteins, lipids)  Organelles Cells Microorganisms Bacteria Phytoplankton neurotoxins Viruses	UV resonance Raman	Microspectroscopy
Cells DNA in chromosomes Pigment in granulocytes & lymphocytes  <u>Major applications</u> Immunoassay Protein identification	Visible	MEDIUM/HIGH  Bioassay reader Cytometry
Tissue: Alzheimers, artery, blood (analytes, LGRL), bone, brain, breast, cervix, colorectal, esophagus, GI, glucose, larynx, prostate, skin, thyroid, tumors  <u>Major applications:</u> Cancer diagnosis Cardiology disease Medical microbiology Tumor detection	NIR	HIGH  Portable instrument with fiber probe vs FT-Raman vs. immunofluorescence staining (tissues)  In vivo Peripheral blood  SERS & endoscopy  Animal/vetnarian Human/clinical

## Biospectroscopy – Examples of Raman spectroscopy - SERS

- In vivo tumor targeting and detection
- Biocompatible and non-toxic nano-particles (pegylated gold) for SERS
- Target tumor biomarkers, for example, epidermal growth factor receptors (EGFR) on human cancer cells and xenograft tumor models


# Biospectroscopy – Examples of Raman spectroscopy - SERS



In vivo tumor targeting and spectroscopic detection and surface-enhanced Raman nanoparticle tags  
*Nature Biotechnology* 26, 83-90 (2008)

## Biospectroscopy – Promise of Enhanced Raman (SERS)

- Subtle changes within biomolecules, such as drug interactions, tissue healing, cosmetics, disease diagnosis
- Intercellular SERS localization and interaction. Identification of drug binding to cells for Drug-DNA and cellular interaction analysis
- Investigation of microorganisms in single cells; yeast cell classifications, single bacterium
- Oxygenation measurements of blood and tissue
- Molecular level cancer detection (cervical, lung, throat and
- Immunoassays using SERS and Raman readers





*Your Photonics Partner*

SPECTROMETERS | LASERS | TOTAL SOLUTIONS

## Contact Information

Michael Kayat, Ph.D.  
Vice President  
B&W Tek, Inc  
19 Shea Way  
Newark, DE 19713  
USA

+1 302 368 7824

[mikek@bwtek.com](mailto:mikek@bwtek.com)